NSK

Electric revolution drive unit for dental use

VIVA MATE G5

OPERATION MANUAL

Please read this Operation Manual carefully before use, and file for future reference.

Powerful Partners®

OM-E0669E

MADE IN JAPAN

CE 0197

Thank you for purchasing the VIVA MATE G5.

Please read this Operation Manual carefully before use to become familiar with operation instructions and care & maintenance. Keep this Operation Manual for future reference.

User and Intended Use : Only qualified Dental Professionals may use this product.

Intended use : A drive unit used in conjunction with devices for cutting and grinding teeth, dentures and crowns and for detaching and removing oral prosthesis and foreign matters. (multiple purposes apply when different attachments are fitted)

- This product is intended as dental equipment. Do not use it for other purposes.

Prohibition : Do not use it in a room with a risk of explosion or in the vicinity of inflammable substances.

Classification of equipment

- Type of protection against electric shock :
 - Class II equipment
- Degree of protection against electric shock :
 - Type BF applied part
- Method of sterilization or disinfection recommended by the manufacture :
 - See 7. Sterilization
- Degree of protection against ingress of water as detailed in the current edition of IEC 60529 :
 - Foot Control: IPX1 (Protection against vertical water drops.)
- Degree of safety of application in the presence of a inflammable anesthetic mixture with air or with oxygen or nitrous oxide :
 - This equipment is not suitable for use in air/ inflammable anesthetic gas or oxygen or nitrogen monoxide (laughing gas) / inflammable anesthetic gas.
- Mode of operation :
 - Intermittent operation

- Contents -

Precautions for handling and operation.....	2
1. Package Contents	5
2. Component Names.....	6
3. Installation and Check before Treatment	9
4. Operation	11
5. Changing Batteries (U1068052).....	13
6. Cleaning.....	15
7. Sterilization.....	15
8. Maintenance.....	15
9. Periodical Maintenance Checks.....	17
10. Troubleshooting.....	17
11. Option Parts List	19
12. Warranty.....	19
13. Disposing Product.....	19
14. Symbols	20
15. Specifications	20
16. EMC Information (Electromagnetic Compatibility Information)	21

Precautions for handling and operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
 DANGER	Hazard that could result in personal death or serious injury if the safety instructions are not correctly followed.
 WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.
 CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

DANGER

- Do not attempt to disassemble the product nor tamper with the mechanism except as recommend by NSK in this Operation Manual, as it may result in injury, electric shock or fire.
- If you should notice a fluid leak from the dry cell batteries, deformation or discoloration of the motor exterior, stop using immediately and contact your Authorized NSK Dealer. Fluid leakage accident, electric shock or fire may result.
- If the battery fluid leaks and gets into your eyes, immediately flush your eyes with water and consult with a medical specialist as it may cause vision loss or blindness.
- Do not throw the battery pack into the fire or heat it. Doing so may result in heat or smoke contamination, explosion, and ignition as insulating material may melt, gas discharge valve and safety mechanism may be damaged, or electrolytic solution may catch fire.
- Do not pierce the battery pack with a needle, hit it with a hammer, or stamp it. Doing so may result in heat or smoke contamination, explosion, and ignition as the battery pack may rupture, deform, or get a short circuited inside.
- If the battery leaks or smells, take it away from fire. Leaked electrolytic solution may catch fire and result in smoke contamination, rupture and fire.

WARNING

- Do not insert or pull out the AC power cord with wet hands as there is a risk of electric shock.
- Do not spill water or other fluid, such as antiseptic solution as this may cause a short circuit resulting in a fire or an electric shock.
- Do not use this product in a room where there is a risk of explosion or in the vicinity of inflammable materials. Do not use on a patient who has been administered laughing gas and in the vicinity of such a patient.
- If the battery fluid leaks and adheres to your skin or clothing, immediately flush it out with water as it may cause damage to your skin.
- Do not pull out or insert the motor cord while the power is ON. If the motor cord is pulled out or plugged in while the motor is rotating, you may be injured by the rotating motor.

⚠ CAUTION


- This product may be used only by professionally trained medical personnel, such as Doctors, Dentists and Dental Hygienists in dental clinics, other medical premises including hospitals, and home visits.
- When operating the product always consider the safety of the patient.
- For irrigation, always use water. If a fluid other than water is used, it may cause clogging.
- Always use the genuine AC adaptor. A non-genuine AC adaptor may result in a failure.
- When the increasing contra angle handpiece is used, always irrigate the water. If irrigation is not performed, it may cause heating.
- Make sure that the motor rotation has stopped before replacing the bur or attaching or removing the handpiece. If the motor is rotating, you may get injured by the bur, or the connecting portion of the handpiece may be damaged, resulting in an injury or a failure.
- Do not direct the handpiece or its lighting directly at the eye of the patient or the operator. Doing so may damage eyes.
- Do not short-circuit the connector terminals of the battery pack. Do not carry or store the battery pack with metal necklaces, hairpins and other metal things. The battery pack may be short-circuited and an excessive electric current may cause heat or smoke contamination, rupture or fire. It may also cause the heating of metal objects, such as wires, necklaces and hair pins.
- Do not allow any impact on to the product. Do not drop the product. Personal injury or damage to the product may result.
- Clean, lubricate and sterilize the handpiece immediately after use. If left with blood, etc. inside the product, the blood may become clotted and cause a failure. Follow the relevant Operation Manual for sterilizing method and maintenance.
- Do not lubricate the motor. Doing so will cause heating and a failure.
- Keep the product in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from dust, salt and sulphur.
- Do not sterilize the control unit, the motor and its cord, the AC adaptor and the handpiece stand.
- If saline or grime has adhered to the Control Unit, turn the power OFF first, then wipe with a damp moist cloth and finally wipe clean thoroughly with a dry cloth.
- To clean the product never uses any solvent such as benzine or thinner.
- Do not disconnect the motor cord from the micromotor.
- When used in a room with electromagnetic interference wave, the operation may be affected. If there is a device emitting such waves, do not use this product in its vicinity. If an ultrasonic generator or an electric scalpel is used nearby, turn OFF the power of this product.
- The use of ACCESSORIES such as cables, with the exception of cables sold by the manufacturer of this product as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of this product.
- VIVA MATE G5 should not be used adjacent to, or stacked with, other equipment. If adjacent or stacked use is necessary, this product should be observed to verify normal operation in the configuration in which it will be used.
- Check for vibration, noise and overheating outside the patient's mouth before use. If any abnormalities are found stop using immediately and contact your dealer.
- Do not use or leave the product in a high-temperature environment such as under strong direct sunlight, in a car under a blazing sun, or near the fire or a stove as it may cause overheating or fire due to a failure of the internal circuit.
- Continuous use of this product over a long period of time may overheat the motor and the handpiece and result in an accident.
- When used beyond the maximum limit of the continuous ON time, the motor may become hot. In such a case, stop using it and wait till it cools down before resuming.

- This product is an intermittent operation device. After three minutes of a continuous operation, leave it OFF for ten minutes.

NOTICE

- Turn off the power after use. If the product is stored for a long time, power off the unit and remove AC adaptor from the product.
- The AC plug is the means to cut off commercial power supply, so make sure that the AC plug can be pulled out from the power outlet without delay in an emergency. In addition, do not place anything within approximately 10 cm of the AC plug.
- Always pull AC power cord by holding adaptor body. Do not pull the AC power cord by holding and pulling its cord. Doing so may break wire in the cord.
- For the handling of the handpiece, see the Operation Manual of the handpiece.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- When a chemical, solvent or disinfectant solution adheres to the control unit or motor cord, immediately wipe it off. If it is left adhered, discoloration or deformation of the product may result.
- If the product is not used for a long period check it is functioning correctly before using on a patient.
- Generally rechargeable batteries constantly discharge little by little even if it is not used. Therefore, we suggest you charge the battery immediately before use.
- This product switches itself off when the voltage reduction in the battery pack is detected. When switched on again after being left for some time, the voltage reduction is sometimes not detected. This is due to battery characteristics and not a fault. The voltage reduction and the remaining capacity of the battery pack do not necessarily correspond to each other. It simply provides an indication.
- No special training is required for this device.
- The end user shall be responsible for any judgment that relates to the application of this product to a patient.
- The maximum motor rotation speed is $30,000\text{min}^{-1}\pm 10\%$.
- Applied part for the patient and the operator is the motor.
- This product is rated Medical Electrical equipment. EMC (Electromagnetic compatibility) is described in the documentation included.
- Installation and use of this product requires special precautions regarding EMC according to the EMC information.
- The product during operation could interfere with the computers, LAN cables in vicinity of use or could cause noise in radio receivers nearby.
- If fluid leakage occurs, thoroughly wipe the fluid adhered to battery mounting parts before inserting new batteries.
- Always use batteries of the same type and manufacturer and replace all three batteries at the same time. Using batteries of different types or mixing new and old or exhausted batteries may cause battery fluid leak or battery rupture.
- When inserting a dry cell battery, do not forcibly press the contact spring in the battery box against the end of the battery's negative terminal as it may cause damage to the battery coating, resulting in a short circuit or battery fluid leakage.
- Do not use rechargeable batteries such as a nickel metal hydride battery and nickel-cadmium battery.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.


1 Package Contents


	Part Name	Quantity
1	Control Unit	1
2	AC Cord	1
3	AC Adaptor	1
4	Motor (Including Motor Cord)	1
5	Handpiece Stand	1
6	O-Ring (Black)	3
7	O-Ring (Blue)	1
	Foot Control (Option)	-

2 Component Names


* Control Unit

Rear Surface

Front Surface

◆ Operation Panel

(1) Power Key

For turning the power ON or OFF.

(2) Rotation Speed Selection Key

For setting the rotation speed.

(3) Rotation Direction Key


For selecting the rotation direction.

(4) ON/OFF Key

For switching ON (rotation)/OFF (stopping rotation) the motor rotation.

◆ Display

Speed Display


- Speed Display
When the motor is not rotating, the set maximum rotation speed is displayed.
When the motor is rotating, the current rotation speed is displayed.
 - When the lighting level is set, the set level is displayed. When there is an error, the corresponding error code is displayed.
 - In an auto-cruise operation, the dot on the bottom right corner flashes.
- Power Display
The status of power is displayed.
 - ➔ AC-IN : Power is being supplied from the AC adaptor.
 - 🔋 CHARGE : The internal battery is being charged.
 - 🔋 Remaining Battery : The battery level is getting low. Charge the battery.
- Rotation Direction Display
Lamp OFF denotes forward rotation. Lamp ON denotes reverse rotation.

CAUTION

- Do not place a heavy object or make an impact (dropping object) on the display of the operation panel, as the display may fail.

* Foot Control (Option)

For turning the motor ON (rotation)/OFF (stop). The rotation speed changes in accordance with how strongly the Foot Control is depressed.

CAUTION

- If you try to turn the power ON while the Foot Control is depressed, the motor will not rotate. If you take your foot off the Foot Control and then press it down again, the motor starts rotating.

3 Installation and Check before Treatment

3-1 To connect the Motor

Align '▲' on the motor cord plug with '▲' on the upper part of the motor socket, and insert the motor cord plug firmly till it clicks. (Fig. 1)

To remove, hold the plug and pull it out straight without twisting.


Fig. 1

3-2 To connect the handpiece (Option)

3-2-1 To connect the handpiece (Option)

When connecting an Optic handpiece, insert the handpiece to the motor insert portion and turn till the locating pin clicks.

For a Non Optic handpiece, insert the handpiece to the motor insert portion. (Fig. 2)


Fig. 2

3-2-2 To remove the handpiece (Option)

Hold the motor and the handpiece separately and disconnect with care.

CAUTION

- When using a straight handpiece, make sure to attach a bur and lock it before starting the rotation.

3-3 To connect irrigation tube (Option)

By combining it with VIVA support 2 (Option), a treatment with irrigation will be possible.

CAUTION

- Prior to use, check the compatibility of the irrigation nozzle of the motor and the irrigation tube. If not connected properly, water may leak. Make sure that the irrigation tube and the irrigation nozzle are securely connected.
- Make sure that the motor and the cord are securely tightened. If not securely tightened, the water may leak.
- Before irrigation, make sure that the handpiece used is suitable for internal irrigation.

3-3-1 Compatibility check of the irrigation nozzle

Use the VIVA support 2 irrigation tube set. If other tubes are used, the tube may break or water may leak from the joining area.

3-3-2 Connection to the irrigation nozzle

Insert the irrigation tube all the way into the irrigation nozzle at the back of the motor.


Fig. 3

CAUTION

- If saline or functional water is used, it may result in corrosion, blockage and damage. Always use water for irrigation.

3-3-3 To check the irrigation (Option)

- 1) Connecting the VIVA support 2.
 - 2) Connecting the Irrigation Tube (Option) to Irrigation Nozzle of the motor.
- * For irrigation method, refer to the Operation Manual of the VIVA support 2.

3-4 To connect Foot Control (Option)

Make sure that the arrow mark on the connection plug of the Foot Control is facing up when the connection plug is inserted to the connection socket of the control unit. (Fig. 4)


Fig. 4

3-5 To connect AC adaptor

Correctly align then insert the AC Power Cord into the AC Power Cord Connector Socket at the back of the Control Unit. (Fig.5)



Fig. 5

4 Operation

4-1 Charging

Charge this product using an AC power source. Prior to use, check the charge status of the battery.

The battery status can be checked by the Remaining Battery display. The lamp flashes when the voltage is low. If this happens, please charge immediately.

- 1) Insert the AC adaptor securely into the power source.
- 2) The AC-IN lamp lights.
- 3) The CHARGE lamp lights and charging starts.
(Battery voltage is measured and, only when it is necessary, charging starts.)
- 4) When the CHARGE lamp goes off, the charging is complete.

CAUTION

- Connect to a commercial power outlet only. If connected to any other power source, there is a risk of failure.
- Before the start of charging, battery voltage is measured automatically. Therefore, it may take several seconds before charging starts. This is not a fault.
- The AC-IN lamp indicates that power is supplied. It remains illuminated after charging is complete.
- Generally it takes approximately 2 hours to charge, though the charging time depends on the use and age of the battery and ambient temperature. An older battery pack tends to have a significantly shorter charging time and usage time.

4-2 Power source

To turn the power ON/OFF, press the Power Key for one second or longer. When the power is ON, the maximum rotation speed currently selected is displayed in the Speed Display.

4-3 Setting the rotation speed

Rotation speed can be adjusted by using ▲ or ▼ of the Rotation Speed Selection Key.

NOTICE

- Rotation speed can be changed when the motor is rotating.
- By holding down Rotation Speed Selection Key, the set rotation speed changes continuously.

4-4 Setting the rotation direction

Every time the Rotation Direction Key is pressed, the rotation direction changes. If the reverse direction has been selected, the Rotation Direction Key lamp is lighted.

- NOTICE**
- Rotation direction can be changed when the motor is rotating.
 - When the Rotation Direction Key is pressed, after a short pause, the rotation in the set direction starts softly.

4-5 Setting the luminance of the lighting

- 1) Turn On the power.
- 2) Press the ON/OFF Key for 2 seconds or longer.
- 3) The luminance setting screen is displayed in the Speed Display.
- 4) Select the luminance using the Rotation Speed Selection Key (▲▼).
 - 1 → Hi mode (bright)
 - 2 → Lo mode (normal)
 - OF → No lighting
- 5) Press the ON/OFF Key to save this setting.

- NOTICE**
- Luminance cannot be set using the Foot Control.

4-6 Operation

When the ON/OFF Key is pressed, the motor rotation starts. When the Key is pressed again, the motor rotation stops (Alternate operation). When Foot Control (Option) operation is used rotation continues as long as the Foot Control is depressed. The rotation speed can be adjusted by depressing the foot control.

4-7 Auto-cruise function (Foot Control operation only)

While the motor is rotating and the rotation speed has been set by the Foot Control, press the ON/OFF key on the control unit to start auto-cruise mode. Auto-cruise mode is indicated by the flashing of the dot on the right bottom corner in the Speed Display. In auto-cruise mode, the rotation and the set speed are maintained even if you take your foot away from the Foot Control.

To cancel auto-cruise mode, either press the ON/OFF key on the control unit briefly or depress the Foot Control again.

- When used beyond the maximum limit of the continuous ON time, the motor may become hot. In such a case, stop using it and wait till it cools down before resuming.
- If the rotation has been started by the ON/OFF Key, auto-cruise function cannot be set up.

4-8 Auto Power Off Function

After 10 minutes of non-activity, power is switched off automatically to save the battery and power consumption. (Auto Power Off function)

To resume operation, turn on the power again.

- NOTICE**
- While the motor handpiece is rotating, the power does not go off even if the keys are not operated.

4-9 At the end of the treatment

When treatment is completed, place the motor handpiece in the holder and press the Power Key for one second or longer to turn OFF the power.

NOTICE • Even if the power is turned off using the Power Key, the set condition is saved, so you do not need to reset again when the power is turned ON next time.

4-10 Adjustment of illumination delay time (0 to 9 seconds: it can be set in multiples of one second.)

Illumination delay time after the rotation has stopped can be adjusted.

- 1) While the power is OFF, press the Power Key while at the same time pressing the ON/OFF key.
- 2) Illumination delay time is displayed on the first digit of the Speed Display.
- 3) Select the suitable illumination delay time (seconds) using the Rotation Speed Selection Key (▲▼).
- 4) Press the ON/OFF Key to save this setting.

5 Changing Batteries (U1068052)

The motor handpiece uses rechargeable batteries. They can be recharged 300-500 times, depending on the use conditions of the motor handpiece. If the operating time becomes shorter or the rotation speed becomes slower, although the "MEMORY EFFECT" described in "NOTICE" is not applicable, the batteries may be at the end of their life expectancy. In such a case, ask your Authorized NSK Dealer to replace the batteries or replace them by yourself with new ones. (Refer to the 5-1 "Changing the batteries".) When replacing batteries by yourself, be sure to observe the following "CAUTIONS ON CHANGING BATTERIES". Note that NSK shall not be liable for any malfunction or failure resulting from you changing the batteries yourself and not following the "CAUTIONS ON CHANGING BATTERIES".

CAUTIONS ON CHANGING BATTERIES

- Do not open any part other than the battery cover.
- Use only batteries as specified by NSK.
- Do not work with wet hands. Failure to do so may result in rust formation on battery terminals or moisture intrusion inside, and could cause failure of the product.
- Do not insert the AC adaptor while the battery is removed from the control unit.
- Should the battery fluid leak and get into your eyes, immediately wash thoroughly with clean water and seek medical attention.

5-1 Changing the batteries

- 1) Turn the power OFF and remove the AC adaptor.
- 2) Remove the screw securing the battery cover on the reverse side of the control unit using a Phillips screwdriver.
- 3) Remove the battery cover by sliding it in the direction of the arrow. (Fig. 6)


Fig. 6

- 4) By pulling out the tape attached to the battery pack, take out the existing battery pack. (Fig. 7)


Fig. 7

- 5) To remove the connector, press the clasp underneath the connector to release it. (Fig.8)

- 6) Insert a new battery securely until its clasp portion locks.
- 7) Insert the battery pack into the control unit, taking care that the cord is not caught.
- 8) Secure the battery cover with the screw.
- 9) Confirm that the power can be turned ON.

Charge the battery before use.


Fig. 8

**CAUTION**

- Do not misplace the screw after removing it.
- Do not throw the battery into fire. The battery may explode and cause an accident.

NOTICE

- The used lithium ion batteries are recyclable, but their disposal may sometimes not be permitted by the country. Return them to your Dealer.

6 Cleaning

Clean the surface of the control unit and motor with a soft brush (do not use a metal brush) and wipe with cotton moistened with rubbing alcohol before and after use. To clean the handpiece, etc., see the enclosed Operation Manual.

**CAUTION**

- Do not lubricate the motor.
- To clean the products never use any solvent such as benzine or thinner.

7 Sterilization

This product is not suitable for sterilization.

*For the sterilization of handpieces, etc., check enclosed Operation Manuals.

8 Maintenance

8-1 Replacing O-Ring

If the O-Ring becomes worn, such as difficult to connect the handpiece or air or water might leak, replace the O-Ring.

Remove the O-Ring at the Motor Insert with a pointed tool, and locate the new O-Rings into the groove.

O-Ring (Blue) : Order Code 0313084070

O-Ring (Black) : Order Code 0312074080



Fig. 9

**CAUTION**

- There are 4 O-Rings on the Motor Insert. The blue O-Ring is thinner than other 3 O-Rings. Make sure to order it.
- Do not use Autoclave Plug without the mounting O-Ring. It may cause a malfunction. If O-Ring had damaged, replace it immediately.
- What can happen if the O-Ring becomes worn
 - Water may leak
 - Water does not come out
 - Vibration
 - Difficult to connect/disconnect the handpiece

8-2 Protection circuit

An electronic circuit breaker automatically functions to protect the motor and the Control Unit if the motor is ever overloaded. Power supply to the motor will automatically be terminated and an Error code will be displayed on the Control Unit.

8-3 Error Code

When the motor handpiece has stopped due to failure, excessive load, disconnection or incorrect use, the control unit condition is automatically programmed to identify the cause of abnormality, and an error code is displayed on the display panel. In such a case, turn the power OFF first and then ON again to see if the same error code is displayed. If the error code is displayed again, refer to the table below for checks and remedy.

Item	Error Code	Error	Cause	Check and Remedy
When motor handpiece is rotating	E.1.	Overcurrent	Motor handpiece is locked.	Remove load.
	E.2.	Overvoltage	Circuit failure	Contact your Authorized NSK Dealer and request repair.
	E.4.	Motor overheating	A high load has been continuously applied on the motor handpiece for a long period of time.	Leave the device for 5 minutes or so till the error code is automatically cancelled. You can then resume using it.
When charging	E.d.	Charging current error	Failure of the battery pack or circuit	Contact your Authorized NSK Dealer and request repair.
	E.E	Outside the scope of operating temperature	Outside the operating temperature range of the battery pack or the failure of the thermistor inside the battery.	Use it within the scope of operating environment (operating conditions). If the error code is still displayed, contact your Authorized NSK Dealer and request repair.

*While an error code is displayed, the CHARGE display flashes.

9 Periodical Maintenance Checks

Perform periodical checks based on the check list below every three months. If abnormality is detected in any of the check items, contact your Authorized NSK Dealer.

As daily Maintenance Checks, follow the Operation Manual and confirm there is no abnormality on Control Unit.

Check item	Description
Rotation	Rotate the motor and check for abnormalities such as vibration, noise and heat.
Rotation speed display	Confirm that the rotation speed is displayed in the Speed Display as usual.
Power display	Insert the AC adaptor and confirm that the AC-IN lamp lights. Also, confirm the CHARGE lamp illuminates when charging.
Rotation direction display	Press the Rotation Direction Key and confirm that Rotation Direction Display lamp lights.
Operation of each button	Press each button and confirm that it functions normally.
Irrigation	Rotate the motor when a handpiece is attached and check for abnormalities such as the rate of irrigation and water leak.

10 Troubleshooting

When a problem is detected, check the following again before requesting a repair. If none of these is applicable or if the trouble is not remedied even after an action has been taken, a failure of this product is suspected. Contact your Authorized NSK Dealer.

<Control unit, AC Adaptor>

Trouble	Cause	Remedy
The power does not turn on.	AC adaptor is not connected.	Check the connection.
	Either the plug of the AC adaptor is not inserted to the power outlet or there is no power supply to the power outlet.	Check the connection.
	The remaining battery is low (when driven by the battery).	Either charge the battery pack or use the AC adaptor.
	Internal fuse has blown.	Request repair.
AC IN lamp does not light.	AC adaptor is not connected.	Check the connection.
	Plug of AC connector is not inserted to the power outlet, or there is no power supply to the power outlet.	Check the connection.
	Internal fuse has blown.	Request repair.
The battery does not charge. (CHARGE lamp does not light.)	Battery is fully charged, or battery is nearly fully charged.	No problems. Continue to use it.
	Battery is charging below the minimum temperature of the operating conditions.	Charge it within the scope of operating environment (operating conditions).
	Battery is charging above the maximum temperature of the operating conditions.	
	An error code is displayed.	Refer to Error Code.
Torque is insufficient or the power supply is suddenly cut off.	The remaining capacity of battery pack is small, or battery pack is not connected (Isn't the Battery Level lamp on and flashing?)	Charge the battery.

< Motor, Handpiece >

Trouble	Cause	Remedy
Motor does not rotate.	Disconnection either inside motor or inside the motor cord.	Request repair from your Authorized NSK Dealer.
	Connectors are not connected properly.	Check that each connector is fully inserted.
Motor does not rotate when Foot Control (Option) is depressed.	Foot Control (Option) is not connected.	Connect Foot Control (Option).
Motor does not rotate when power is turned ON.	Handpiece is clogged up.	Clean or replace handpiece.
	Short circuit inside the handpiece or in the motor cord.	Request repair from your Authorized NSK Dealer.
	Power is turned ON while the ON/OFF Key is pressed.	Take your foot away from Foot Control (Option), then operate with ON/OFF Key.
	Short circuit inside ON/OFF Key.	Request repair from your Authorized NSK Dealer.
	Power is turned ON while Foot Control (Option) is depressed.	Check Foot Control (Option).
	Short circuit inside Foot Control (Option).	Disconnect Foot Control (Option) and operate with the ON/OFF Key, or contact your dealer and request repair from your Authorized NSK Dealer.
LED lamps do not light.	Connectors are not connected properly.	Check that each connector is fully inserted.
	Disconnection either inside motor or inside the motor.	Request repair from your Authorized NSK Dealer.
	LED is broken.	Request repair from your Authorized NSK Dealer.
	Luminance setting is OFF.	Select either 1 or 2 in the luminance setting. (See 4-5 Setting the luminance of the lighting.)
Abnormal vibration, heat, noise, abnormal smell, faulty rotation of motor	There is damage either inside the motor or in the motor cord.	Request repair from your Authorized NSK Dealer.
Water leaks from motor, or no water comes out in irrigation.	The connectors are not connected properly.	Check that each connector is fully inserted.
	O-Ring has deteriorated.	Replace the O-Ring. (See 8-1 Replacing O-Ring).
	VIVA support 2 (Option) does not operate.	Check the Operation Manual of VIVA support 2 (Option).
	The rate of irrigation from VIVA support 2 (Option) is not sufficient.	VIVA MATE G5 motor is an internal irrigation motor. So, when used with VIVA support 2, the irrigation rate may be insufficient. In such a case, please consult your Authorized NSK Dealer.

Water leaks from motor, or no water comes out in irrigation.	Irrigation tube (Option) is not adequate.	Use a compatible irrigation tube.
	Motor inside or motor code inside is broken.	Request repair from your Authorized NSK Dealer.
Rotation cannot be stopped by Foot Control (Option).	Rotation has been activated by the ON/OFF Key.	Use the ON/OFF Key to stop rotation.

*Servicing, such as checking the inside of the product, is performed by service personnel according to the service manual, so contact your Authorized NSK Dealer for servicing as needed.

11 Option Parts List

Model name	Model	Order Code
Battery Pack		U1068052
O-Ring (Blue)		0313084070
O-Ring (Black)		0312074080
Foot Control (FC-76M)		Z1082002

12 Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

13 Disposing Product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

The used batteries are recyclable, but their disposal may sometimes not be permitted by the respective country.

14 Symbols

  Consult operation instructions.

 Class II Equipment

 Type BF applied part.

 CAUTION

  Conforms to CE European Directive of "Medical equipment directive 93/42/EEC."

  Mandatory conformance mark on many products placed on the market in the European Economic Area

 Follow the waste of electric and electronic equipment (WEEE) Directive (2002/96/EC) for product and accessory disposal.

 Manufacturer.

 Authorized representative in the European community

 For indoor use only.

 Efficiency Level

15 Specifications

< Control Unit >

Model	NE280
Rated Input	DC24V 2.5A
Rated Output	15W
Charging time	About 2 H
Dimensions	W72×D180×H36mm
Weight	370g

< AC adaptor , AC Cord >

Model	TR60M24
Rated Input	AC100 - 240V
Dimensions	W58×D132×H31mm
Weight	410g (Include AC cord)

< Motor >

Model	VMG30EM
Rated Input	DC12V 1.3A
Dimensions	Ø22×L68.6mm
Weight	148g (Include Motor cord)
Rotation	3,000 - 30,000min ⁻¹
Torque	3N · cm (30,000 min ⁻¹ while rotating)
Optic	White LED

	Temperature	Humidity	Atmospheric Pressure
Use environment	0 - 40°C (32 - 104°F) (No Condensation)	30 - 75% RH	/
Store and transportation environment	-10 - 50°C (14 - 122°F)	10 - 85% RH	

*No liquid freezing allowed

**When operated outside the above range, there is a risk of failure.

16 EMC Information (Electromagnetic Compatibility Information)

Guidance and manufacturer's declaration - Electromagnetic Emissions		
VIVA MATE G5 is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR11/ EN55011	Group 1	VIVA MATE G5 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11/EN55011	Class B	VIVA MATE G5 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions EN/IEC61000-3-2	Not Applicable	
Voltage fluctuations/flicker emissions EN/IEC 61000-3-3	Not Applicable	

Guidance and manufacturer's declaration - Electromagnetic Immunity			
VIVA MATE G5 is intended for use in the electromagnetic environment specified below. The customer or the user of VIVA MATE G5 should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) EN/IEC61000-4-2	±(2,4,6)kV contact ±(2,4,8)kV air	±(2,4,6)kV contact ±(2,4,8)kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EN/IEC61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN/IEC61000-4-5	±1kV line(s) to line(s) ±2kV line(s) to earth	±1kV line(s) to line(s) ±2kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN/IEC61000-4-11	<5% Ut (>95% dip in Ut) for 0.5 cycles 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 sec	<5% Ut (>95% dip in Ut) for 0.5 cycles 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of VIVA MATE G5 requires continued operation during power mains interruptions, it is recommended that VIVA MATE G5 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz)magnetic field EN/IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: 'Ut' is the AC mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - Electromagnetic Immunity

VIVA MATE G5 is intended for use in the electromagnetic environment specified below. The customer or the user of VIVA MATE G5 should assure that it is used in such an environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF EN/IEC61000-4-6	3Vrms 150kHz to 80MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.11 \sqrt{P}$ $d = 1.00 \sqrt{P}$ 80MHz to 800MHz $d = 2.00 \sqrt{P}$ 800MHz to 2.5GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, and (d) is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey ⁽⁸⁾ should be less than the compliance level in each frequency range ⁽⁸⁾ . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF EN/IEC61000-4-3	3V/m 80MHz to 2.5GHz	3 V/m	

NOTE1: At 80MHz and 800MHz, the higher frequency range applies.

NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the product is used exceeds the applicable RF compliance level stated above, the product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product.

b: Over the 150kHz to 80MHz frequency range, the field strength should be less than 3V/m.

c: Externally displayed symbol of a device or component of device comprising an RF transmitter or emitting RF magnetic field energy for diagnosis or treatment. 

Cable and Accessories	Max length	Complies with	
Motor Cord	1.5m	RF emission CISPR11	Class B/ Group 1
Foot Control Cord	2.0m		
AC Adaptor	1.8m	Electrostatics discharge(ESD) Surge immunity Power frequency magnetic field Conducted RF Radiated RF	EN/IEC61000-4-2 EN/IEC61000-4-5 EN/IEC61000-4-8 EN/IEC61000-4-6 EN/IEC61000-4-3

Recommended separation distances between portable and mobile RF communications equipment and VIVA MATE G5

VIVA MATE G5 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of VIVA MATE G5 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and VIVA MATE G5 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz - 80MHz $d = 1.2 \sqrt{P}$	80MHz - 800MHz $d = 1.2 \sqrt{P}$	800MHz - 2.5GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance 'd' in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where 'P' is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1: At 80MHz and 800MHz, the higher frequency range applies.

NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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NSK Europe GmbH
NSK France SAS
NSK United Kingdom Ltd
NSK Asia Pte Ltd