

THIS IS A NEXT-GENERATION STANDARD MACHINE.

NSK has sought best-in-class performance, operability and durability with ULTIMATE 500.

A thoroughgoing reduction of costs has realized excellent cost performance.

With special consideration given to safety, it has cleared the world's rigorous safety standards.


It has been designed eco-friendlily, minimizing the power consumption.

The advanced design includes a stylish form.

NSK aims at the top without compromising on every point.



ULTIMATE 500



UM500LGC
Light Grey / Compact set
NE95-LG
UM50C
FC-40

ULTIMATE 500 Set & Spare Parts

■ ULTIMATE 500 SET

Product	Set	Code No.	
UM500LGT Light Grey / Torque set	NE95-LG	Y140 551 (120V)	
	UM50T FC-40	Y140 552 (230V)	
UM500LGC Light Grey / Compact set	NE95-LG	Y140 554 (120V)	
	UM50C FC-40	Y140 555 (230V)	
UM500SBT Translucent Blue / Torque set	NE95-SB	Y140 557 (120V)	
	UM50T FC-40	Y140 558 (230V)	
UM500SBC Translucent Blue / Compact set	NE95-SB	Y140 560 (120V)	
	UM50C FC-40	Y140 561 (230V)	

■ SPARE PARTS

Product	Code No.	
Control Unit (with Battery Cord and Handpiece Stand)	NE95-LG Light Grey	U385 001 (120V) U386 001 (230V)
	NE95-SB Translucent Blue	U385 002 (120V) U386 002 (230V)
Motor Handpieces	UM50T Torque type	Y140 562
	UM50C Compact type	Y140 563
Motor (without Cord)	UM50TM Torque type	E113
	UM50CM Compact type	E187
Handpieces	UHR50T Torque type	H165
	UHR50C Compact type	H189
Motor Cord	Torque type	E113 051
	Compact type	E187 051
Foot Pedal	FC-40	Z056
Handpiece Stand (without Spanner)		Z095 201
Chuck	Ø 2.35 CHH2.35	K003 180
	Ø 1.6 CHH1.6	K003 183
	Ø 3.0 CHH3.0	K003 179



NSK, an environmentally friendly company, has received recognition of its conformity to the Environment Management System "ISO 14001", the International Standard of quality assurance "ISO 9001" and approval of using CE Markings (MDD) from the German Qualification Organization TÜVRheinland.

- Always read the instruction manual before operating any new product.
- Specifications may vary from country to country and are subject to change without notice.

NSK NAKANISHI INC.

HEAD OFFICE & FACTORY : 340 Kamihinata, Kanuma-shi, Tochigi-ken 322-8666, Japan
TEL: 0289-64-3380 FAX: 0289-62-5636
TOKYO OFFICE : 3F Sakai Bldg., Ueno 3-19-4, Taito-ku, Tokyo 110-0005, Japan
TEL: 03-3835-2892 FAX: 03-3835-2856

<http://www.nsk-nakanishi.co.jp> e-mail: box@nsk-nakanishi.co.jp